


CIL NMR Solvent Data Chart

More Solvents, More Sizes,
More Solutions

Solvent	¹ H Chemical Shift (ppm from TMS) (multiplicity)	JHD (Hz)	Carbon-13 Chemical Shift (ppm from TMS) (multiplicity)	JCD (Hz)	¹ H Chemical Shift of HOD (ppm from TMS)	Density at 20°C	Melting point (°C)	Boiling point (°C)	Dielectric Constant	Molecular Weight
Acetic Acid-d ₄	11.65 (1) 2.04 (5)	2.2	178.99 (1) 20.0 (7)	20	11.5	1.12	16.7	118	6.1	64.08
Acetone-d ₆	2.05 (5)	2.2	206.68 (1) 29.92 (7)	0.9 19.4	2.8	0.87	-94	56.5	20.7	64.12
Acetonitrile-d ₃	1.94 (5)	2.5	118.69 (1) 1.39 (7)	21	2.1	0.84	-45	81.6	37.5	44.07
Benzene-d ₆	7.16 (1)		128.39 (3)	24.3	0.4	0.95	5.5	80.1	2.3	84.15
Chloroform-d	7.24 (1)		77.23 (3)	32.0	1.5	1.50	-63.5	61-62	4.8	120.38
Cyclohexane-d ₁₂	1.38 (1)		26.43 (5)	19	0.8	0.89	6.47	80.7	2.0	96.24
Deuterium Oxide	4.80 (DSS) 4.81 (TSP)		NA	NA	4.8	1.11	3.81	101.42	78.5	20.03
N, N-Dimethyl-formamide-d ₇	8.03 (1) 2.92 (5) 2.75 (5)	1.9 1.9	163.15 (3) 34.89 (7) 29.76 (7)	29.4 21.0 21.1	3.5	1.03	-61	153	36.7	80.14
Dimethyl Sulfoxide-d ₆	2.50 (5)	1.9	39.51 (7)	21.0	3.3	1.19	18.45	189	46.7	84.17
1,4-Dioxane-d ₈	3.53 (m)		66.66 (5)	21.9	2.4	1.13	11.8	101.1	2.2	96.16
Ethanol-d ₆	5.19 (1) 3.56 (1) 1.11 (m)		56.96 (5) 17.31 (7)	22 19	5.3	0.89	-114.1	78.5	24.5	52.11
Methanol-d ₄	4.78 (1) 3.31 (5)	1.7	49.15 (7)	21.4	4.9	0.89	-97.8	64.7	32.7	36.07
Methylene Chloride-d ₂	5.32 (3)	1.1	54.00 (5)	27.2	1.5	1.35	-95	39.75	8.9	86.95
Pyridine-d ₅	8.74 (1) 7.58 (1) 7.22 (1)		150.35 (3) 135.91 (3) 123.87 (3)	27.5 24.5 25	5	1.05	-42	115-116	12.4	84.13
1,1,2,2-Tetrachloroethane-d ₄	6.0		73.78 (3)			1.62	-44	147	8.20	169.86
Tetrahydrofuran-d ₈	3.58 (1) 1.73 (1)		67.57 (5) 25.37 (5)	22.2 20.2	2.4-2.5	0.99	-108.5	66	7.6	80.16
Toluene-d ₈	7.09 (m) 7.00 (1) 6.98 (5) 2.09 (5)	2.3	137.86 (1) 129.24 (3) 128.33 (3) 125.49 (3) 20.4 (7)	23 24 24 19	0.4	0.94	-95	110.6	2.4	100.19
Trifluoroacetic Acid-d	11.50 (1)		164.2 (4) 116.6 (4)		11.5	1.41	-15.4	72.4		115.03
Trifluoroethanol-d ₃	5.02 (1) 3.88 (4x3)	2(9)	126.3 (4) 61.5 (4x5)	22	5	1.41	-43.3	75		103.06

S Budavari, M.J. O'Neil, A. Smith, P.E. Heckelman, *The Merck Index, an Encyclopedia of Chemicals, Drugs, and Biologicals - Eleventh Edition*, Merck Co., Inc. Rahway, NJ, 1989.



Cambridge Isotope Laboratories, Inc.

50 Frontage Road, Andover, MA 01810-5413 USA

PH: 800.322.1174 (N. AMERICA) PH: 978.749.8000

FAX: 978.749.2768 WEB: www.isotope.com

EMAIL: cilsales@isotope.com (N. AMERICA)

EMAIL: intl@isotope.com (INTERNATIONAL)


CIL

The ¹H spectra of the residual protons and ¹³C spectra were obtained on a Varian Gemini 200 spectrometer at 295°K. The NMR solvents used to acquire these spectra contain a maximum of 0.05% and 1.0% TMS (v/v) respectively. Since deuterium has a spin of 1, triplets arising from coupling to deuterium have the intensity ratio of 1:1:1. 'm' denotes a broad peak with some fine structures. It should be noted that chemical shifts can be dependent on solvent, concentration and temperature.

Approximate values only, may vary with pH, concentration and temperature.

Melting and boiling points are those of the corresponding unlabeled compound (except for D₂O). These temperature limits can be used as a guide to determine the useful liquid range of the solvents. Information gathered from the Merck Index - Eleventh Edition.

NMR Solvent Storage and Handling Information



Acetic Acid-d₄	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Acetone-d₆	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Acetonitrile-d₃	Store at room temperature away from light and moisture. Stable for one year after receipt of order if stored under above conditions (unopened). After one year, the compound should be re-checked for chemical purity before use.
Benzene-d₆	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Chloroform-d	Store refrigerated (-5° to 5°C) away from light and moisture. Stable for six months after receipt of order if stored under above conditions (unopened). After six months, the compound should be re-checked for chemical purity before use.
Cyclohexane-d₁₂	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Deuterium Oxide	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
N,N-Dimethylformamide-d₇	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Dimethyl Sulfoxide-d₆	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
1,4-Dioxane-d₈ (p-Dioxane)	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Ethanol-d₆	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Methanol-d₄	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Methylene Chloride-d₂	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Pyridine-d₅	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
1,1,2,2-Tetrachloroethane-d₂	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Tetrahydrofuran-d₈	Store refrigerated (-5° to 5°C) away from light and moisture. Stable for six months after receipt of order if stored under above conditions (unopened). After six months, the compound should be re-checked for chemical purity before use.
Toluene-d₈	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
Trifluoroacetic Acid-d	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.
2,2,2-Trifluoroethanol-d₃	Store at room temperature away from light and moisture. Stable if stored under recommended conditions.



Cambridge Isotope Laboratories, Inc.

50 Frontage Road, Andover, MA 01810-5413 USA
 PH: 800.322.1174 (N. AMERICA) PH: 978.749.8000
 FAX: 978.749.2768 WEB: www.isotope.com
 EMAIL: cilsales@isotope.com (N. AMERICA)
 EMAIL: intlsales@isotope.com (INTERNATIONAL)

Please note that some packaging sizes of some solvents may require special handling not given above. The bottle or ampule packaging should be checked for further instructions.